

WWW.SAFE-DIETS.COM

DEFINITION

GELDIET Transport is a complete nutritional solution that can be used to both diet and hydrate animals using a single product.



Photo indicative only

PRODUCT OBJECTIVE

GELDIET Transport is composed of 59.2 % water and of proteins, slow carbohydrates (starch), lipids, vitamins and minerals, to provide animals with their main nutrients during transport or a stressful situation such as a post-operative period.

GELDIET Transport can be distributed provisionally (2 to 7 days) to animals in a stressful situation or a recovery phase as a complement to granulated diets.

DIRECTIONS FOR USE

Distribution period:

- 1 day before the transport phase for acclimatization, then throughout transport.
- 72 hours for animals weakened by protocol conditions; replenishment varies depending on cage type (particularly cage ventilation).

GELDIET Transport can be distributed

- ✓ Placed directly on the cage floor, in the open package,
- ✓ Removed from the packaging and placed directly in the cage dieting dish, or on the cage floor.

Quantity consumed per day: 10-15 g / day for mice, 40-50 g / day for rats. Daily consumption depends on species and strain.

Distribution method: ad libitum or rationing depending on experimental protocols. Standard diet (e.g. A04 or A03 granules) can be kept within animals' reach.

WWW.SAFE-DIETS.COM

Storage:

The product in its closed, original packaging can be stored for 6 months in a dark area at room temperature. Can be stored at 4°C.

After opening, the product can be kept for 5 days maximum.

ABOUT THIS PRODUCT

Solid paste with 59.2% water content, packaged in individual, protectively sealed tubs. Irradiated at 10 kiloGrays.



Photo indicative only

- ✓ 100 g tub, sold in packs of 60. The tubs are grouped in quantities of 30 in a hermetically sealed envelope.
- ✓ 200 g tub, sold in packs of 36. The tubs are grouped in quantities of 18 in a hermetically sealed envelope.

NUTRITIONAL COMPOSITION

Calorie content (kcal/kg) 1,300

✓ Water	59.2%
✓ Proteins	8.5%
✓ Lipids	1.6%
✓ Carbohydrates	21.5%
✓ Fibre	1.5%

INGREDIENTS

Water, sugars, casein, maltodextrin, cellulose, lard, vitamin, minerals, hydrocolloids, preservatives.