

**WWW.SAFE-DIETS.COM**

## DEFINITION

GELDIET Classic is a complete nutritional solution that can be used to both diet and hydrate animals using a single product.



Photo indicative only

## PRODUCT OBJECTIVE

GELDIET Classic is composed of 73% water and 25% A03 or A04-type standard diet. This Geldiet provides a wide range of nutrients: proteins, slow carbohydrates (starch), lipids, vitamins and minerals. It can be used to provide animals with their main nutrients during transport or a stressful situation such as a post-operative period.

GELDIET Classic can be distributed provisionally (2 to 7 days) to animals in a stressful situation or a recovery phase as a complement to granulated diets.

## DIRECTIONS FOR USE

### Distribution period:

- 1 day before the transport phase for acclimatization, then throughout transport.
- 72 hours for animals weakened by protocol conditions; replenishment varies depending on cage type (particularly cage ventilation).

### GELDIET Classic can be distributed :

- ✓ Placed directly on the cage floor, in the open package
- ✓ Removed from the packaging and placed directly in the cage dieting dish, or on the cage floor.

**Quantity consumed per day:** 10-15 g/day for mice, 40-50 g/day for rats. Daily consumption depends on species and strain.

**Distribution method :** ad libitum or rationing depending on experimental protocols. Standard diet (e.g. A04 or A03 granules) can be kept within animals' reach.

**[WWW.SAFE-DIETS.COM](http://WWW.SAFE-DIETS.COM)**

## Storage:

The product in its closed, original packaging can be stored for 6 months in a dark area at room temperature. Can be stored at 4°C.

After opening, the product can be kept for 5 days maximum.

## ABOUT THIS PRODUCT

Solid paste with 73% water content, packaged in individual, protectively sealed tubs. Irradiated at 10 kilograys.

- ✓ 100 g tub, sold in packs of 60. The tubs are grouped in quantities of 30 in a hermetically sealed envelope.
- ✓ 200 g tub, sold in packs of 36. The tubs are grouped in quantities of 18 in a hermetically sealed envelope.

## NUTRITIONAL COMPOSITION

Calorie content (kcal/kg) 1,300

- |                 |     |
|-----------------|-----|
| ✓ Water         | 73% |
| ✓ Standard diet | 25% |

## INGREDIENTS

Water, standard diet (A03 or A04), hydrocolloids, preservatives.